

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse COL28A1 Polyclonal Antibody

多克隆抗体

产品货号: PA04499

第三版

描述:COL28A1 (Collagen Type XXVIII Alpha 1 Chain) is a Protein Coding gene. Among its related pathways are Collagen chain trimerization and Degradation of the extracellular matrix. GO annotations related to this gene include serine-type endopeptidase inhibitor activity. An important paralog of this gene is COL6A2. COL28A1 belongs to a class of collagens containing von Willebrand factor (VWF; MIM 613160) type A (VWFA) domains (Veit et al., 2006 [PubMed 16330543]).

抗原:Synthesized peptide derived from Collagen alpha-1(XXVIII)Polyclonal Antibody

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: