

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat IMP3 Polyclonal Antibody

多克隆抗体

产品货号: PA04527

第三版

**描述:**This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain. IMP3 (IMP3, U3 Small Nucleolar Ribonucleoprotein) is a Protein Coding gene. Diseases associated with IMP3 include Enchondroma. Among its related pathways are MAPK-Erk Pathway and rRNA processing in the nucleus and cytosol. GO annotations related to this gene include poly(A) RNA binding and structural constituent of ribosome.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**