

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MIB1 Polyclonal Antibody

多克隆抗体

产品货号: PA04530

第三版

描述:MIB1 (Mindbomb E3 Ubiquitin Protein Ligase 1) is a Protein Coding gene. Diseases associated with MIB1 include Left Ventricular Noncompaction 7 and Left Ventricular Noncompaction. Among its related pathways are Notch signaling pathway (KEGG) and Signaling by NOTCH1. GO annotations related to this gene include ligase activity and ubiquitin-protein transferase activity. An important paralog of this gene is MIB2. This gene encodes a protein containing multiple ankyrin repeats and RING finger domains that functions as an E3 ubiquitin ligase. The encoded protein positively regulates Notch signaling by ubiquitinating the Notch receptors, thereby facilitating their endocytosis. This protein may also promote the ubiquitination and degradation of death-associated protein kinase 1 (DAPK1).

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed at low level. Expressed at higher level in spinal cord, ovary, whole brain, and all specific brain regions examined.