

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MYOZ2 Polyclonal Antibody

多克隆抗体

产品货号: PA04547

第三版

描述:MYOZ2 (Myozenin 2) is a Protein Coding gene. Diseases associated with MYOZ2 include Cardiomyopathy, Hypertrophic, 16 and Myoz2-Related Familial Hypertrophic Cardiomyopathy. GO annotations related to this gene include actin binding and telethonin binding. An important paralog of this gene is MYOZ1. The protein encoded by this gene belongs to a family of sarcomeric proteins that bind to calcineurin, a phosphatase involved in calcium-dependent signal transduction in diverse cell types. These family members tether calcineurin to alpha-actinin at the z-line of the sarcomere of cardiac and skeletal muscle cells, and thus they are important for calcineurin signaling. Mutations in this gene cause cardiomyopathy familial hypertrophic type 16, a hereditary heart disorder.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed specifically in heart and skeletal muscle.