

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat FAM48A Polyclonal Antibody

多克隆抗体

产品货号: PA04548

---

第三版

**描述:**SUPT20H (SPT20 Homolog, SAGA Complex Component) is a Protein Coding gene. Among its related pathways are Autophagy - animal and Chromatin organization. GO annotations related to this gene include transcription cofactor activity.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Highly expressed in testis, moderately in brain and pituitary gland. Expressed in several fetal tissues, including lung, brain, thymus and kidney. Expression is down-regulated in malignant prostate tissues.