
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CHID1 Polyclonal Antibody

多克隆抗体

产品货号: PA04553

第三版

描述:CHID1 (Chitinase Domain Containing 1) is a Protein Coding gene. Among its related pathways are Response to elevated platelet cytosolic Ca²⁺. GO annotations related to this gene include hydrolase activity, hydrolyzing O-glycosyl compounds and chitinase activity.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in cells of monocytic, T, B and epithelial origin.