

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD54 Polyclonal Antibody

多克隆抗体

产品货号: PA04589

第三版

**描述:** ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). During leukocyte trans-endothelial migration, ICAM1 engagement promotes the assembly of endothelial apical cups through ARHGEF26/SGEF and RHOG activation. In case of rhinovirus infection acts as a cellular receptor for the virus.

**抗原:** Synthesized peptide derived from human ICAM-1 around the non-phosphorylation site of Tyr512.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**