

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MTIF3 Polyclonal Antibody

多克隆抗体

产品货号: PA04593

第三版

描述: This gene encodes a translation initiation factor that is involved in mitochondrial protein synthesis. Polymorphism in this gene is associated with the onset of Parkinson's disease. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 5.

抗原: Synthesized peptide derived from the Internal region of human IF3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: