

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IFITM3 Polyclonal Antibody

多克隆抗体

产品货号: PA04595

第三版

描述:The protein encoded by this gene is an interferon-induced membrane protein that helps confer immunity to influenza A H1N1 virus, West Nile virus, and dengue virus. Two transcript variants, only one of them protein-coding, have been found for this gene. Another variant encoding an N-terminally truncated isoform has been reported, but the full-length nature of this variant has not been determined.

抗原:Synthesized peptide derived from the N-terminal region of human IFITM3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: