

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IGFBP1 Polyclonal Antibody

多克隆抗体

产品货号: PA04597

第三版

描述:This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors.

抗原:Synthesized peptide derived from the Internal region of human IGFBP1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: