

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IgG1 Polyclonal Antibody

多克隆抗体

产品货号: PA04599

第三版

描述:IGHG1 (Immunoglobulin Heavy Constant Gamma 1 (G1m Marker)) is a Protein Coding gene. Diseases associated with IGHG1 include Chronic Lymphocytic Leukemia and Heavy Chain Deposition Disease. Among its related pathways are Role of phospholipids in phagocytosis and Innate Immune System. GO annotations related to this gene include antigen binding. An important paralog of this gene is IGHG3.

抗原:Synthesized peptide derived from the Internal region of human IgG1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: