

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL32 Polyclonal Antibody

多克隆抗体

产品货号: PA04610

第三版

描述: Cytokine that may play a role in innate and adaptive immune responses. It induces various cytokines such as TNFA/TNF-alpha and IL8. It activates typical cytokine signal pathways of NF-kappa-B and p38 MAPK.

抗原: Synthesized peptide derived from the C-terminal region of human IL-32.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Selectively expressed in lymphocytes. Expression more prominent in immune cells than in non-immune cells.