

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human INSL4 Polyclonal Antibody

多克隆抗体

产品货号: PA04615

第三版

描述:INSL4 encodes the insulin-like 4 protein, a member of the insulin superfamily. INSL4 encodes a precursor that undergoes post-translational cleavage to produce 3 polypeptide chains, A-C, that form tertiary structures composed of either all three chains, or just the A and B chains. Expression of INSL4 products occurs within the early placental cytotrophoblast and syncytiotrophoblast. INSL4 (Insulin Like 4) is a Protein Coding gene. GO annotations related to this gene include receptor binding and insulin-like growth factor receptor binding.

抗原:Synthesized peptide derived from the Internal region of human INSL4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in placenta, uterus and in fetal perichondrium. Expression levels were increased in both early placentas and molar pregnancies and were reduced in choriocarcinoma cells.