

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ITPKC Polyclonal Antibody

多克隆抗体

产品货号: PA04616

第三版

描述: This gene encodes a member of the inositol 1,4,5-trisphosphate [Ins(1,4,5)P(3)] 3-kinase family of enzymes that catalyze the phosphorylation of inositol 1,4,5-trisphosphate to 1,3,4,5-tetrakisphosphate. The encoded protein is localized to the nucleus and cytoplasm and has both nuclear import and nuclear export activity. Single nucleotide polymorphisms in this gene are associated with Kawasaki disease. ITPKC (Inositol-Trisphosphate 3-Kinase C) is a Protein Coding gene. Diseases associated with ITPKC include Kawasaki Disease and Lymph Node Disease. Among its related pathways are Inositol phosphate metabolism (KEGG) and Calcium signaling pathway. GO annotations related to this gene include calmodulin binding and inositol-1,4,5-trisphosphate 3-kinase activity. An important paralog of this gene is ITPKB.

抗原: Synthesized peptide derived from the Internal region of human InsP 3-kinase C

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in pancreas, skeletal muscle, liver, placenta and weakly in kidney and brain.