

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse IP6K3 Polyclonal Antibody

多克隆抗体

产品货号: PA04617

第三版

描述:This gene encodes a protein that belongs to the inositol phosphokinase (IPK) family. This protein is likely responsible for the conversion of inositol hexakisphosphate (InsP6) to diphosphoinositol pentakisphosphate (InsP7/PP-InsP5). It may also convert 1,3,4,5,6-pentakisphosphate (InsP5) to PP-InsP4. Alternative splicing results in multiple transcript variants encoding the same protein. IP6K3 (Inositol Hexakisphosphate Kinase 3) is a Protein Coding gene. Among its related pathways are Farnesoid X Receptor Pathway and Inositol phosphate metabolism (REACTOME). GO annotations related to this gene include inositol-1,4,5-trisphosphate 3-kinase activity and inositol hexakisphosphate 1-kinase activity. An important paralog of this gene is IP6K1.

抗原:Synthesized peptide derived from the Internal region of human InsP6 Kinase 3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in brain.