

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat ICK Polyclonal Antibody

多克隆抗体

产品货号: PA04619

---

第三版

**描述:**ICK may play a key role in the development of multiple organ systems and particularly in cardiac development.

**抗原:**Synthesized peptide derived from human Intestinal Cell Kinase around the non-phosphorylation site of Y159.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in heart, brain, placenta, pancreas, thymus, prostate, testis, ovary, small intestine and colon, with highest levels in placenta and testis. Not detected in spleen. Also expressed in many cancer cell lines.