

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IVL Polyclonal Antibody

多克隆抗体

产品货号: PA04620

第三版

**描述:**Involucrin, a component of the keratinocyte crosslinked envelope, is found in the cytoplasm and crosslinked to membrane proteins by transglutaminase. This gene is mapped to 1q21, among calpactin I light chain, trichohyalin, profilaggrin, loricrin, and calyculin.

**抗原:**Synthesized peptide derived from the C-terminal region of human Involucrin

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Keratinocytes of epidermis and other stratified squamous epithelia.