

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat JDP2 Polyclonal Antibody

多克隆抗体

产品货号: PA04645

第三版

**描述:**JDP2 (Jun Dimerization Protein 2) is a Protein Coding gene. Diseases associated with JDP2 include Keratitis and Skeletal Muscle Cancer. Among its related pathways are Tacrolimus/Cyclosporine Pathway, Pharmacodynamics and ATF-2 transcription factor network. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is ATF3.

**抗原:**Synthesized peptide derived from human JDP2 around the non-phosphorylation site of Thr148.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**