

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse JMY Polyclonal Antibody

多克隆抗体

产品货号: PA04649

第三版

**描述:**JMY (Junction Mediating And Regulatory Protein, P53 Cofactor) is a Protein Coding gene. Diseases associated with JMY include Usher Syndrome, Type 1B. Among its related pathways are Regulation of TP53 Activity and Gene Expression. GO annotations related to this gene include actin binding and transcription coactivator activity. An important paralog of this gene is WHAMM.

**抗原:**Synthesized peptide derived from the C-terminal region of human JMY

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**