



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PPP1R14C Polyclonal Antibody

多克隆抗体

产品货号: PA04663

第三版

描述:PP14C is an inhibitor of PPP1CA. It has over 600-fold higher inhibitory activity when phosphorylated, creating a molecular switch for regulating the phosphorylation status of PPP1CA substrates and smooth muscle contraction. It is found in breast carcinoma.

抗原:Synthesized peptide derived from the Internal region of human KEPI

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in breast cancer.