

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human COL12A1 Polyclonal Antibody

多克隆抗体

产品货号: PA04667

第三版

描述:COL12A1 (Collagen Type XII Alpha 1 Chain) is a Protein Coding gene. Diseases associated with COL12A1 include Bethlem Myopathy 2 and Ullrich Congenital Muscular Dystrophy 2. Among its related pathways are Collagen chain trimerization and Phospholipase-C Pathway. GO annotations related to this gene include extracellular matrix structural constituent conferring tensile strength. An important paralog of this gene is COL14A1. This gene encodes the alpha chain of type XII collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XII collagen is a homotrimer found in association with type I collagen, an association that is thought to modify the interactions between collagen I fibrils and the surrounding matrix. Alternatively spliced transcript variants encoding different isoforms have been identified.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Found in collagen I-containing tissues: both isoform 1 and isoform 2 appear in amnion, chorion, skeletal muscle, small intestine, and in cell culture of dermal fibroblasts, keratinocytes and endothelial cells. Only isoform 2 is found in lung, placenta, kidney and a squamous cell carcinoma cell line. Isoform 1 is also present in the corneal epithelial Bowman's membrane (BM) and the interfibrillar matrix of the corneal stroma, but it is not detected in the limbal BM.