

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat KCNJ2 Polyclonal Antibody

多克隆抗体

产品货号: PA04668

第三版

**描述:** Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein, which has a greater tendency to allow potassium to flow into a cell rather than out of a cell, probably participates in establishing action potential waveform and excitability of neuronal and muscle tissues. Mutations in this gene have been associated with Andersen syndrome, which is characterized by periodic paralysis, cardiac arrhythmias, and dysmorphic features.

**抗原:** Synthesized peptide derived from the Internal region of human KIR2.1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Heart, brain, placenta, lung, skeletal muscle, and kidney. Diffusely distributed throughout the brain.