

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SUV39H2 Polyclonal Antibody

多克隆抗体

产品货号: PA04681

第三版

**描述:**SUV39H2 (Suppressor Of Variegation 3-9 Homolog 2) is a Protein Coding gene. Among its related pathways are Circadian rhythm related genes and Lysine degradation. GO annotations related to this gene include chromatin binding and transcription regulatory region sequence-specific DNA binding. An important paralog of this gene is SUV39H1.

**抗原:**Synthesized peptide derived from the Internal region of human KMT1B

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**