

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TCL1 Polyclonal Antibody

多克隆抗体

产品货号: PA04688

第三版

**描述:**TCL1A (T-Cell Leukemia/Lymphoma 1A) is a Protein Coding gene. Diseases associated with TCL1A include Precursor T-Cell Acute Lymphoblastic Leukemia and T-Cell Leukemia. Among its related pathways are PI3K / Akt Signaling and PI3K-Akt signaling pathway. An important paralog of this gene is MTCP1.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Restricted in the T-cell lineage to immature thymocytes and activated peripheral lymphocytes. Preferentially expressed early in T- and B-lymphocyte differentiation.