

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LMTK2 Polyclonal Antibody

多克隆抗体

产品货号: PA04690

第三版

**描述:**The protein encoded by this gene belongs to the protein kinase superfamily and the protein tyrosine kinase family. It contains N-terminal transmembrane helices and a long C-terminal cytoplasmic tail with serine/threonine/tyrosine kinase activity. This protein interacts with several other proteins, such as Inhibitor-2 (Inh2), protein phosphatase-1 (PP1C), p35, and myosin VI. It phosphorylates other proteins, and is itself also phosphorylated when interacting with cyclin-dependent kinase 5 (cdk5)/p35 complex. This protein involves in nerve growth factor (NGF)-TrkA signalling, and also plays a critical role in endosomal membrane trafficking. Mouse studies suggested an essential role of this protein in spermatogenesis. LMTK2 (Lemur Tyrosine Kinase 2) is a Protein Coding gene. Diseases associated with LMTK2 include Prostate Cancer Susceptibility and Deafness, Autosomal Recessive 37. Among its related pathways are GPCR Pathway and Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is AATK.

**抗原:**Synthesized peptide derived from the Internal region of human KPI-2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Mainly expressed in skeletal muscle, and weakly in brain and pancreas.