
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RPA3 Polyclonal Antibody

多克隆抗体

产品货号: PA04696

第三版

描述:RPA3 (Replication Protein A3) is a Protein Coding gene. Diseases associated with RPA3 include Epidemic Typhus and Ovary Serous Adenocarcinoma. Among its related pathways are DNA Double-Strand Break Repair and Telomere C-strand (Lagging Strand) Synthesis. GO annotations related to this gene include single-stranded DNA binding and damaged DNA binding.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: