

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GADD34 Polyclonal Antibody

多克隆抗体

产品货号: PA04700

第三版

**描述:**PPP1R15A (Protein Phosphatase 1 Regulatory Subunit 15A) is a Protein Coding gene. Diseases associated with PPP1R15A include Malignant Pleural Mesothelioma. Among its related pathways are Transcriptional activity of SMAD2/SMAD3-SMAD4 heterotrimer and Activation of cAMP-Dependent PKA. GO annotations related to this gene include protein kinase binding. This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The induction of this gene by ionizing radiation occurs in certain cell lines regardless of p53 status, and its protein response is correlated with apoptosis following ionizing radiation.

**抗原:** Synthetic peptide from human protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:**