

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat KCNS2 Polyclonal Antibody

多克隆抗体

产品货号: PA04709

第三版

描述:KCNS2 (Potassium Voltage-Gated Channel Modifier Subfamily S Member 2) is a Protein Coding gene. Among its related pathways are Potassium Channels and Transmission across Chemical Synapses. GO annotations related to this gene include ion channel activity and delayed rectifier potassium channel activity. An important paralog of this gene is KCNS3.

抗原:Synthesized peptide derived from the Internal region of human KV9.2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: