

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SH2D2A Polyclonal Antibody

多克隆抗体

产品货号: PA04714

第三版

描述: This gene encodes an adaptor protein thought to function in T-cell signal transduction. A related protein in mouse is responsible for the activation of lymphocyte-specific protein-tyrosine kinase and functions in downstream signaling. Alternative splicing results in multiple transcript variants. SH2D2A (SH2 Domain Containing 2A) is a Protein Coding gene. Among its related pathways are Development VEGF signaling via VEGFR2 - generic cascades and Immune response Fc epsilon RI pathway. GO annotations related to this gene include SH3 domain binding and SH3/SH2 adaptor activity. An important paralog of this gene is SH2D7.

抗原: Synthesized peptide derived from the Internal region of human Lad

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expression limited to tissues of the immune system and, in particular, activated T-cells. Expressed in peripheral blood leukocytes, thymus and spleen. Much lower expression or undetectable, in brain, placenta, skeletal muscle, prostate, testis, ovary, small intestine, and colon. Expressed at low levels in unstimulated T-cells, but not expressed in normal resting or activated B-cells. According to PubMed:10692392, expression is not restricted to activated T-cells, but strongly expressed in blood cell lineages, the endothelium and other cell and tissue types, such as heart, lung, and liver.