

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DDT Polyclonal Antibody

多克隆抗体

产品货号: PA04763

第三版

**描述:**DDT (D-Dopachrome Tautomerase) is a Protein Coding gene. Diseases associated with DDT include Disuse Amblyopia. GO annotations related to this gene include dopachrome isomerase activity and D-dopachrome decarboxylase activity. An important paralog of this gene is DDTL.D-dopachrome tautomerase converts D-dopachrome into 5,6-dihydroxyindole. The DDT gene is related to the migration inhibitory factor (MIF) in terms of sequence, enzyme activity, and gene structure. DDT and MIF are closely linked on chromosome 22.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**