

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IFNAR1 Polyclonal Antibody

多克隆抗体

产品货号: PA04764

第三版

描述:IFNAR1 (Interferon Alpha And Beta Receptor Subunit 1) is a Protein Coding gene. Diseases associated with IFNAR1 include Septicemic Plague and Hepatitis C. Among its related pathways are Type I Interferon Signaling Pathways and Immune response IFN alpha/beta signaling pathway. GO annotations related to this gene include interferon receptor activity and type I interferon receptor activity. An important paralog of this gene is ENSG00000249624. The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:IFN receptors are present in all tissues and even on the surface of most IFN-resistant cells. Isoform 1, isoform 2 and isoform 3 are expressed in the IFN-alpha sensitive myeloma cell line U266S. Isoform 2 and isoform 3 are expressed in the IFN-alpha resistant myeloma cell line U266R. Isoform 1 is not expressed in IFN-alpha resistant myeloma cell line U266R.