

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ITGA1 Polyclonal Antibody

多克隆抗体

产品货号: PA04771

第三版

描述:ITGA1 (Integrin Subunit Alpha 1) is a Protein Coding gene. Diseases associated with ITGA1 include Phimosis and Colon Carcinoma In Situ. Among its related pathways are GPCR Pathway and MAPK-Erk Pathway. GO annotations related to this gene include receptor binding and collagen binding. An important paralog of this gene is ITGA2. This gene encodes the alpha 1 subunit of integrin receptors. This protein heterodimerizes with the beta 1 subunit to form a cell-surface receptor for collagen and laminin. The heterodimeric receptor is involved in cell-cell adhesion and may play a role in inflammation and fibrosis. The alpha 1 subunit contains an inserted (I) von Willebrand factor type I domain which is thought to be involved in collagen binding.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: