

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LRG1 Polyclonal Antibody

多克隆抗体

产品货号: PA04778

第三版

描述:The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002 [PubMed 12223515]).

抗原:Synthesized peptide derived from the C-terminal region of human LRG1.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Plasma.