

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human INSL3 Polyclonal Antibody

多克隆抗体

产品货号: PA04782

第三版

描述:INSL3 (Insulin Like 3) is a Protein Coding gene. Diseases associated with INSL3 include Cryptorchidism and Testicular Disease. Among its related pathways are Peptide ligand-binding receptors and Signaling by GPCR. GO annotations related to this gene include receptor binding and protease binding. This gene encodes a member of the insulin-like hormone superfamily. The encoded protein is mainly produced in gonadal tissues. Studies of the mouse counterpart suggest that this gene may be involved in the development of urogenital tract and female fertility. This protein may also act as a hormone to regulate growth and differentiation of gubernaculum, and thus mediating intra-abdominal testicular descent. Mutations in this gene may lead to cryptorchidism. Alternate splicing results in multiple transcript variants.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in prenatal and postnatal Leydig cells. Found as well in the corpus luteum, trophoblast, fetal membranes and breast.