

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GH Polyclonal Antibody

多克隆抗体

产品货号: PA04807

第三版

**描述:**GH1 (Growth Hormone 1) is a Protein Coding gene. Diseases associated with GH1 include Growth Hormone Deficiency, Isolated, Type Ii and Kowarski Syndrome. Among its related pathways are GPCR Pathway and Phospholipase-C Pathway. GO annotations related to this gene include growth factor activity and growth hormone receptor binding. An important paralog of this gene is CSH1. GH2 (Growth Hormone 2) is a Protein Coding gene. Diseases associated with GH2 include Choriocarcinoma and Prostate Carcinoma In Situ. Among its related pathways are GPCR Pathway and Phospholipase-C Pathway. GO annotations related to this gene include hormone activity. An important paralog of this gene is GH1. The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. The gene, along with four other related genes, is located at the growth hormone locus on chromosome 17 where they are interspersed in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. The five genes share a remarkably high degree of sequence identity. Alternative splicing generates additional isoforms of each of the five growth hormones, leading to further diversity and potential for specialization.

**抗原:** Synthetic peptide from human protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**