

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey MARCH3 Polyclonal Antibody

多克隆抗体

产品货号: PA04828

第三版

描述: This gene encodes a member of the membrane-associated RING-CH (MARCH) family. The encoded protein is an E3 ubiquitin-protein ligase that may be involved in regulation of the endosomal transport pathway. MARCH3 (Membrane Associated Ring-CH-Type Finger 3) is a Protein Coding gene. GO annotations related to this gene include ligase activity. An important paralog of this gene is MARCH2.

抗原: Synthesized peptide derived from the Internal region of human March3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: