

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MARK1/2/3/4 Polyclonal Antibody

多克隆抗体

产品货号: PA04831

第三版

描述: May play a role in cytoskeletal stability.

抗原: Synthesized peptide derived from the Internal region of human MARK1/2/3/4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: MARK1 is Highly expressed in heart, skeletal muscle, brain, fetal brain and fetal kidney. MARK2 is High levels of expression in heart, brain, skeletal muscle and pancreas, lower levels observed in lung, liver and kidney.1 Publication MARK3 and MARK4 are Ubiquitous.