

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat/Monkey MARK1/2/3/4 Polyclonal Antibody

多克隆抗体

产品货号: PA04832

---

第三版

**描述:** May play a role in cytoskeletal stability.

**抗原:** Synthesized peptide derived from human MARK1/2/3/4 around the non-phosphorylation site of Thr215.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** MARK1 is Highly expressed in heart, skeletal muscle, brain, fetal brain and fetal kidney. MARK2 is High levels of expression in heart, brain, skeletal muscle and pancreas, lower levels observed in lung, liver and kidney. MARK3 and MARK4 are Ubiquitous.