
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MAST2 Polyclonal Antibody

多克隆抗体

产品货号: PA04841

第三版

描述:MAST2 (Microtubule Associated Serine/Threonine Kinase 2) is a Protein Coding gene. Diseases associated with MAST2 include Breast Mucinous Carcinoma. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is MAST4.

抗原:Synthesized peptide derived from the Internal region of human MAST205

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Abundant in the testis.