

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MAST3 Polyclonal Antibody

多克隆抗体

产品货号: PA04843

第三版

描述:MAST3 (Microtubule Associated Serine/Threonine Kinase 3) is a Protein Coding gene. Diseases associated with MAST3 include Inflammatory Bowel Disease 6. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is MAST2.

抗原:Synthesized peptide derived from the C-terminal region of human MAST3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: