

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ST14 Polyclonal Antibody

多克隆抗体

产品货号: PA04851

第三版

**描述:**Degrades extracellular matrix. Proposed to play a role in breast cancer invasion and metastasis. Exhibits trypsin-like activity as defined by cleavage of synthetic substrates with Arg or Lys as the P1 site.

**抗原:**Synthesized peptide derived from the N-terminal region of human Matrilysin

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**