

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat M-Cadherin Polyclonal Antibody

多克隆抗体

产品货号: PA04867

第三版

描述: This gene is a member of the cadherin superfamily of genes, encoding calcium-dependent intercellular adhesion glycoproteins. Cadherins consist of an extracellular domain containing 5 cadherin domains, a transmembrane region, and a conserved cytoplasmic domain. Transcripts from this particular cadherin are expressed in myoblasts and upregulated in myotubule-forming cells. The protein is thought to be essential for the control of morphogenetic processes, specifically myogenesis, and may provide a trigger for terminal muscle cell differentiation. CDH15 (Cadherin 15) is a Protein Coding gene. Diseases associated with CDH15 include Autosomal Dominant Non-Syndromic Intellectual Disability and 16Q24.3 Microdeletion Syndrome. Among its related pathways are Cell junction organization and Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CDH4.

抗原: Synthesized peptide derived from the N-terminal region of human M-cadherin

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in the brain and cerebellum.