

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MDFI Polyclonal Antibody

多克隆抗体

产品货号: PA04875

第三版

**描述:** This protein is a transcription factor that negatively regulates other myogenic family proteins. Studies of the mouse homolog, I-mf, show that it interferes with myogenic factor function by masking nuclear localization signals and preventing DNA binding. Knockout mouse studies show defects in the formation of vertebrae and ribs that also involve cartilage formation in these structures. MDFI (MyoD Family Inhibitor) is a Protein Coding gene. GO annotations related to this gene include transcription factor binding. An important paralog of this gene is MDFIC.

**抗原:** Synthesized peptide derived from the Internal region of human MDFI

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:**