

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MDM1 Polyclonal Antibody

多克隆抗体

产品货号: PA04880

第三版

**描述:** This gene encodes a nuclear protein similar to the mouse double minute 1 protein. The mouse gene is located in double minute (DM) chromatin particles, is amplified in the mouse transformed 3T3 cell line, and the encoded protein is able to bind to p53. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. MDM1 (Mdm1 Nuclear Protein) is a Protein Coding gene. Diseases associated with MDM1 include Retinitis Pigmentosa 3.

**抗原:** Synthesized peptide derived from the C-terminal region of human MDM1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**