

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HDGF Polyclonal Antibody

多克隆抗体

产品货号: PA04893

第三版

描述:HDGF,also named as HMG1L2,belongs to the HDGF family. It is a heparin-binding protein,with mitogenic activity for fibroblasts. HDGF acts as a transcriptional repressor. It plays important roles in the development and progression of cancers. HDGF may have a synergistic effect with VEGF during cancer development and Progression.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.