

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MED8 Polyclonal Antibody

多克隆抗体

产品货号: PA04894

第三版

描述: This gene encodes a protein component of the mediator complex, which aids in transcriptional activation through interaction with RNA polymerase II and gene-specific transcription factors. The encoded protein may also function in ubiquitin ligation and protein degradation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. MED8 (Mediator Complex Subunit 8) is a Protein Coding gene. Among its related pathways are Influenza A and Regulation of lipid metabolism by Peroxisome proliferator-activated receptor alpha (PPARalpha). GO annotations related to this gene include RNA polymerase II transcription cofactor activity.

抗原: Synthesized peptide derived from the Internal region of human Med8

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: