

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GSTA1 Polyclonal Antibody

多克隆抗体

产品货号: PA04908

第三版

描述: GSTA1 (Glutathione S-Transferase Alpha 1) is a Protein Coding gene. Diseases associated with GSTA1 include Ovarian Endodermal Sinus Tumor and Ovarian Primitive Germ Cell Tumor. Among its related pathways are Glutathione metabolism and Platinum drug resistance. GO annotations related to this gene include glutathione transferase activity. An important paralog of this gene is GSTA2. This gene encodes a member of a family of enzymes that function to add glutathione to target electrophilic compounds, including carcinogens, therapeutic drugs, environmental toxins, and products of oxidative stress. This action is an important step in detoxification of these compounds. This subfamily of enzymes has a particular role in protecting cells from reactive oxygen species and the products of peroxidation. Polymorphisms in this gene influence the ability of individuals to metabolize different drugs. This gene is located in a cluster of similar genes and pseudogenes on chromosome 6. Alternative splicing results in multiple transcript variants.

抗原: Synthetic peptide from human protein

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Liver.