

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GSC Polyclonal Antibody

多克隆抗体

产品货号: PA04924

第三版

描述:GSC (Goosecoid Homeobox) is a Protein Coding gene. Diseases associated with GSC include Short Stature, Auditory Canal Atresia, Mandibular Hypoplasia, Skeletal Abnormalities and Auditory Neuropathy. Among its related pathways are Regulation of retinoblastoma protein and Embryonic and Induced Pluripotent Stem Cell Differentiation Pathways and Lineage-specific Markers. GO annotations related to this gene include sequence-specific DNA binding and RNA polymerase II repressing transcription factor binding. An important paralog of this gene is GSC2. This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: