

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat COL14A1 Polyclonal Antibody

多克隆抗体

产品货号: PA04933

第三版

**描述:**COL14A1 (Collagen Type XIV Alpha 1 Chain) is a Protein Coding gene. Diseases associated with COL14A1 include Keratoderma, Palmoplantar, Punctate Type Ia and Punctate Palmoplantar Keratoderma. Among its related pathways are Collagen chain trimerization and Phospholipase-C Pathway. GO annotations related to this gene include poly(A) RNA binding and extracellular matrix structural constituent. An important paralog of this gene is COL12A1. This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**