

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MERTK Polyclonal Antibody

多克隆抗体

产品货号: PA04938

第三版

**描述:** This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP).

**抗原:** Synthesized peptide derived from the N-terminal region of human MerTK

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Not expressed in normal B- and T-lymphocytes but is expressed in numerous neoplastic B- and T-cell lines.